Puritan Regan-Lowe Semi-Solid Transport Medium

Intended Use

Puritan Regan-Lowe Transport Medium is a semi-solid medium intended for transport of specimens containing *Bordetella* species.

Summary and Explanation

Pertussis (whooping cough) is a respiratory infection that can be severe in young infants caused by fastidious bacteria such as *Bordetella pertussis* and *Bordetella parapertussis*. Early diagnosis and treatment are essential to limit the progression of the disease and minimize its transmission. Regan-Lowe Transport Medium is a selective medium used for transport of nasopharyngeal specimens to the laboratory for qualitative procedures. The base medium contains charcoal to neutralize toxic fatty acids in the specimen. It is supplemented with cephalexin to inhibit bacteria indigenous to the nasopharynx. Defibrinated horse blood is added to the medium to support the growth of *Bordetella* species.^{1,2}

Summary

Beef extract and gelatin provide amino acids and nitrogen to support bacteria growth. Sodium chloride helps to maintain osmotic balance. Defibrinated horse blood and niacin provide nutrients for growth of *Bordetella* species. Cephalexin is added to inhibit normal flora of the nasopharynx.³

Approximate formula

Beef Extract.....	5.0 g	Niacin.....	0.005g
Pancreatic digest of gelatin.....	5.0 g	Defibrinated Horse Blood.....	100 mL
Starch.....	5.0 g	Cephalexin.....	40 mg
Charcoal.....	2.0 g	Agar.....	6.0 g
Sodium chloride.....	2.5 g	Demineralized Water.....	1000 mL

PH 7.4 ± 0.2 @ 25°C

Procedures

1. Collect specimen from the posterior nasopharynx with a swab or aspirate.
2. Place swab or specimen into Regan-Lowe transport medium and transport to the laboratory for culture.
3. Once in the laboratory inoculate onto an enrichment selective Regan-Lowe Agar plate.
4. Incubate aerobically at 35°-36°C at least for 7 days.^{3,4} To avoid drying agar plates should be placed in a moist chamber with a damp filter paper.
5. Check plates daily for growth. *Bordetella* colonies appear small, smooth, transparent and glistening on Regan-Lowe Agar plates.³

Quality control

All batches of Puritan Regan-Lowe Transport Medium are tested prior to release for pH and further evaluated for their ability to promote growth of the following organism:

CONTROL

Bordetella pertussis ATCC 9797

INCUBATION

Aerobic, 48 h @ Room Temperature

RESULTS

Good Recovery

Limitation

For specimen collection and transport only. Once in the laboratory, specimen should be streaked onto agar plate for culture and subsequent biochemical tests. Specimens collected after administration of antibiotics may have lower rate of bacteria recovery.

References


1. Regan, J., F. Lowe. 1977. Enrichment Medium for the Isolation of *Bordetella*. J. Clin. Microbiol. 6(3): 303-309.
2. Hoppe, J.E., M. Schlagenhauf. 1989. Comparison of Three Kinds of Blood and Two Incubation Atmospheres for Cultivation of *Bordetella pertussis* on Charcoal Agar. J. Clin. Microbiol. 27(9): 2115-2117.
3. Zimbro, M.J., D.A. Power, S.M. Miller, G.E. Wilson, J.A. Johnson. 2009. Difco & BBL Manual of Microbiological Culture Media, 2nd ed. Becton Dickinson and Company. Sparks, MD.
4. Murphy, T., K. Bisgard, G. Sanden. 2000. Diagnosis and Laboratory Methods. <http://www.cdc.gov/pertussis/outbreaks/guide/downloads/chapter-02-amended.pdf>. Accessed January 7, 2016.

207-876-3311 • puritanmedproducts.com

sales@puritanmedproducts.com

Puritan Medical Products Co. LLC
31 School Street, Guilford, Maine 04443-0149 USA
ISO 9001:2008 ISO 13485:2003 CE

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