## Puritan CT Broth-transportmedium

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1. Rivers, C.A., J. R. Schwebke. 2008. Viability of *Trichomonas vaginalis* in Copan Universal Transport Medium and eSwab Transport Medium. *Journal of Clinical Microbiology*. 46(9): 3134-3135.
2. Zimbro M.J., D.A. Power. 2003. Difco & BBL Manual: Manual of Microbiological Culture Media. Becton, Dickinson, and Company. Sparks, MD.
3. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risk related exposure to biological agents at work. *Official Journal of the European Communities*. L 262/21-45.
4. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, D.W. Warnock. 2011. *Manual of Clinical Microbiology*, 10<sup>th</sup> ed. American Society for Microbiology. Washington, DC.
5. Miller, J.M. 1996. *A guide to specimen management in clinical microbiology*. American Society for Microbiology. Washington, DC.
6. Forbes, B.A., D.F. Sahm, A.S. Weissfeld. 2007. *Diagnostic Microbiology* 12<sup>th</sup> ed. Mosby. St. Louis, MO.
7. Murray, P.R., E.G. Baron, J.H. Jorgensen, M.A. Pfaller, R.H. Tenover, R.M. Tenover. 2003. *Manual of Clinical Microbiology*, 8<sup>th</sup> ed. American Society for Microbiology, Washington, DC.
8. CLSI. *Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard-Third Edition*. CLSI document M22-A3. Wayne, PA. Clinical and Laboratory Standards Institute; 2004.

207-876-3311 • [puritanmedproducts.com](http://puritanmedproducts.com)  
[sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)  
Puritan Medical Products Co. LLC  
31 School Street, Guilford, Maine 04443-0149 USA  
ISO 9001:2008 ISO 13485:2003 CE