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1. Jones, D.E., E.M. Friedl, K.S. Kanarek, J.K. Williams, and D.V. Lim. 1983. Rapid identification of pregnant women heavily colonized with group B streptococci. *J Clin Microbiol.* 18:558-560.
2. Jones, D.E., K.S. Kanarek, D.V. Lim. 1984. Group B streptococcal colonization patterns in mothers and their infants. *J Clin Microbiol.* 20:438-440.
3. Church, D.L., H. Baxter, T. Lloyd, B. Miller, S. Elsayed. 2008. Evaluation of StrepB Carrot Broth versus Lim Broth for detection of Group B streptococcus colonization status of near-term pregnant women. *J Clin Microbiol.* 46(8):2780-2782.
4. Schrag, S, R. Gorwitz, K. Fultz-Butts, A. Schuchat. 2002. Prevention of Perinatal Group B Streptococcal Disease. *Morbidity and Mortality Weekly Report* 51:1-26.
5. Verani, J.R., L. McGee, S.J. Schrag. 2010. Prevention of Perinatal Group B Streptococcal Disease. *Morbidity and Mortality Weekly Report* 59:1-32.
6. Elsayed S., D.B. Gregson, D.L. Church. 2003. Comparison of direct selective versus nonselective agar media plus Lim Broth enrichment for determination of Group B streptococcus colonization status in pregnant women. *Archives of Pathology & Laboratory Medicine.* 127(6): 718-720.
7. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, D.W. Tenover, W. Tenover. 2011. *Manual of Clinical Microbiology*, 10<sup>th</sup> ed. American Society for Microbiology. Washington, DC.
8. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risk related exposure to biological agents at work. *Official Journal of the European Communities.* L 262/21-45.
9. Miller, J.M. 1996. *A guide to specimen management in clinical microbiology.* American Society for Microbiology. Washington, DC.
10. Forbes, B.A., D.F. Sahm, A.S. Weissfeld. 2007. *Diagnostic Microbiology* 12<sup>th</sup> ed. Mosby. St. Louis, MO.
11. Murray, P.R., E.G. Baron, J.H. Jorgensen, M.A. Pfaller, R.H. Tenover, R. Tenover. 2003. *Manual of clinical microbiology*, 8<sup>th</sup> ed. American Society for Microbiology, Washington, D.C.
12. CLSI. *Quality Control for Commercially Prepared Microbiological Culture Media; Approved Standard-Third Edition.* CLSI document M22-A3. Wayne, PA. Clinical and Laboratory Standards Institute; 2004.
13. Zimbro M.J, D.A. Power. 2003. *Difco & BBL Manual: Manual of Microbiological Culture Media.* Becton, Dickinson and Company. Sparks, MD.

207-876-3311 • [puritanmedproducts.com](http://puritanmedproducts.com)  
[sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

Puritan Medical Products Co. LLC  
31 School Street, Guilford, Maine 04443-0149 EE. UU.  
ISO 9001:2008 ISO 13485:2003 ©