## Puritan Regan-Lowe Semi-Solid Transport Medium

### Intended Use

Puritan Regan-Lowe Transport Medium is a semi-solid medium intended for transport of specimens containing *Bordetella* species.

### Summary and Explanation

Pertussis (whooping cough) is a respiratory infection that can be severe in young infants caused by fastidious bacteria such as *Bordetella pertussis* and *Bordetella parapertussis*. Early diagnosis and treatment are essential to limit the progression of the disease and minimize its transmission. Regan-Lowe Transport Medium is a selective medium used for transport of nasopharyngeal specimens to the laboratory for qualitative procedures. The base medium contain charcoal to neutralize toxic fatty acids in the specimen. It is supplemented with cephalixin to inhibit bacteria indigenous to the nasopharynx. Defibrinated horse blood is added to the medium to support the growth of *Bordetella* species.<sup>1,2</sup>

### Summary

Beef extract and gelatin provide amino acids and nitrogen to support bacteria growth. Sodium chloride helps to maintain osmotic balance. Defibrinated horse blood and niacin provide nutrients for growth of *Bordetella* species. Cephalixin is added to inhibit normal flora of the nasopharynx.<sup>3</sup>

### Approximate formula

Beef Extract.....	5.0 g	Niacin.....	0.005g
Pancreatic digest of gelatin.....	5.0 g	Defibrinated Horse Blood.....	100 mL
Starch.....	5.0 g	Cephalixin.....	40 mg
Charcoal.....	2.0 g	Agar.....	6.0 g
Sodium chloride.....	2.5 g	Demineralized Water.....	1000 mL

PH 7.4 ± 0.2 @ 25°C

### Procedures

1. Collect specimen from the posterior nasopharynx with a swab or aspirate.
2. Place swab or specimen into Regan-Lowe transport medium and transport to the laboratory for culture.
3. Once in the laboratory inoculate onto an enrichment selective Regan-Lowe Agar plate.
4. Incubate aerobically at 35°-36°C at least for 7 days.<sup>3,4</sup> To avoid drying agar plates should be placed in a moist chamber with a damp filter paper.
5. Check plates daily for growth. *Bordetella* colonies appear small, smooth, transparent and glistening on Regan-Lowe Agar plates.<sup>3</sup>

### Quality control

All batches of Puritan Regan-Lowe Transport Medium are tested prior to release for pH and further evaluated for their ability to promote growth of the following organism:

#### CONTROL

*Bordetella pertussis* ATCC 9797

#### INCUBATION

Aerobic, 48 h @ Room Temperature

#### RESULTS

Good Recovery

### Limitation


For specimen collection and transport only. Once in the laboratory, specimen should be streaked onto agar plate for culture and subsequent biochemical tests. Specimens collected after administration of antibiotics may have lower rate of bacteria recovery.

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Puritan Medical Products Co. LLC