



ESK - Buffered Peptone Water

ESK® Sampling Kit - Pre-filled with 4 ml or 10 ml Buffered Peptone Water



Description

The Buffered Peptone Water is a pre-enriched non-selective buffered solution. It allows for the repair of damaged cells and facilitate the recovery of target bacteria in samples taken for sanitation monitoring in industrial, food and dairy plants.

Dimensions

4.06" Overall Length of Kit (103.124 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.450 mm)

Handle

- 3.25" Length (82.55 mm)

Qualities

- Easy edge-peel packaging
- Quick-turn leak-resistant cap
- High-quality Puritan polyester tipped swab
- Swab is securely attached to cap
- Shelf-life: 18 months
- Made in USA

Applications

- For environmental surface sampling
- Quality control
- Detection of Microorganisms

Ref #	Packaging
25-83004 PD BPW Sterile	50/1 (4 ml - 50 individually wrapped kits with polyester swab)
25-83004 PDB BPW Sterile	6/50/1 (4 ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)
25-83010 PD BPW Sterile	50/1 (10 ml - 50 individually wrapped kits with polyester swab)
25-83010 PDB BPW Sterile	6/50/1 (10 ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)



Puritan®