



ESK - Butterfield's Solution

ESK Sampling Kit - Pre-filled with 4 ml or 10 ml Butterfield's Solution, 0.1% Peptone and 2% Tween 80



Description

Butterfield's Solution is specified for use in diluting water, dairy products, and food for microbiological methods.

Dimensions

4.06" Overall Length of Kit (103.124 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.450 mm)

Handle

- 3.25" Length (82.55 mm)

Qualities

- Easy edge-peel packaging
- Quick-turn leak-resistant cap
- High-quality Puritan polyester-tipped swab
- Swab is securely attached to cap
- Shelf-life: 18 months
- Made in USA

Applications

- For environmental surface sampling
- Quality control
- Detection of microorganisms

Ref #	Packaging
25-83004 PD BS Sterile	50/1 (4 ml - 50 individually wrapped kits with polyester swab)
25-83004 PDB BS Sterile	6/50/1 (4 ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)
25-83010 PD BS Sterile	50/1 (10 ml - 50 individually wrapped kits with polyester swab)
25-83010 PDB BS Sterile	6/50/1 (10 ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)



Puritan®