



# Flat Stick

Wood, square end



-Choose the shape best for your application

-Good and sturdy, made to work hard

## Dimensions

6" Overall Length (152.4 mm)

Width: 0.75" (19.05 mm)

Thickness: 0.063" (1.6 mm)

## Qualities

- Flat blade
- One end straight cut
- Made in USA

## Applications

- Health and Beauty - waxing
- For mixing and applying materials
- Paints, inks, chemicals
- Stirring stick
- Ideal for crafts
- Other general use

Ref #	Packaging
705 SQ END Non-sterile	10/500 (10 boxes per case, 500 per box)



**Puritan®**