



HydraFlock® Flocked Swab - Dry Transport System

Standard tip, polystyrene handle



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. The patented HydraFlock® flocked swab provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

Tip

- 0.200" Dia. (5.08 mm)
- 0.69" Length (17.526 mm)

Handle

- 0.100" Dia. (2.54 mm) tapered to 0.066" (1.676 mm)
- 6.01" Length (152.654 mm)

Qualities

- Medical grade quality
- Highly absorbent capacity
- Packaged in a reclosable dry transport tube for ease in transport to lab
- Double swab to accommodate multiple sampling sites
- Made in USA

Applications

- For collection and release of biological specimens
- Buccal cell collection
- Specimen collection
- Forensic evidence collection
- Rapid diagnostic collection
- DNA testing

Ref #	Packaging
25-3306-2H BT Sterile	10/50/2 (10 boxes per case of 50 individually wrapped tubes per box, 2 swabs per tube)



Puritan®