

PURITAN MEDICAL PRODUCTS COMPANY, LLC

Safety Data Sheet

Rev. 004 March 10, 2015

1. Product and Company Identification	
Product Name: Puritan UniTranz-RT™ Transport System	
Product Number: UT-117	
Company Identification: Puritan Medical Products Company, LLC P.O. Box 149, 31 School Street Guilford, Maine 04443-0149 U.S.A.	Contact Numbers: Tel: +1 207-876-3311 Fax: +1 207-876-3130
2. Hazards Identification Classification according to Directive 671548/EEC or Directive 1999/45/EC	
Skin contact:	Irritation
Eye contact:	Irritation
3. Composition/information on ingredients	
Composition of Puritan UniTranz-RT™ Transport System consists of the following: Hank's Balance Salt Solution, Bovine Serum Albumin, Gelatin, Sucrose, L-cysteine, , L-glutamic acid, Colistin, Amphotericin B, Vancomycin, Sodium Bicarbonate, and HEPES. All materials listed are considered non-hazardous at ≤1% concentration of final solution.	
4. First-aid measures	
Skin contact:	Immediately wash with water and soap, rinse thoroughly.
Eye contact:	Flush open eye with plenty of water for 15 minutes.
Inhalation:	Consult Doctor.
Swallowing:	Consult Doctor.
5. Fire-fighting measures	
Extinguishing media:	Dry chemicals, CO ₂ , Foam or halogen extinguisher.
Protective Equipment:	Personal Protective Equipment (PPE)
6. Accidental release measures	
Personal precautions:	Personal Protective Equipment (PPE)
Environmental precautions:	N/A
Methods for cleaning up:	N/A
7. Handling and storage	
Handling:	No special measures required
Storage:	Store away from oxidizing agents Store at 2°- 25°C.
8. Exposure controls/personal protection	
Respiratory protection:	N/A
Hand protection:	Disposable Rubber Gloves
Eye protection:	Safety Glasses
Skin and body protection:	Protective work clothing (lab coat)

Safety Data Sheet

9. Physical and chemical properties	
pH:	7.3 ± 0.2
Density:	Not determined
Boiling point, °C:	Not determined
Melting point, °C:	Not determined
Flash point, °C:	Not applicable
Solubility:	Insoluble
10. Stability and reactivity	
Materials and conditions to avoid:	No decomposition if used according to specifications
Hazardous decomposition products:	No dangerous decomposition products known
11. Toxicology information	
Skin	Irritant to skin and mucous membranes
Eye	Irritant effect
Sensitization	No sensitization effect known
12. Ecological Information	
Ecology Effects:	None have been identified.
General Notes:	Not known to be hazardous to water.
13. Disposal considerations	
Recommendation:	Dispose of materials in accordance with federal (40 CFR 261.3), state, and local requirements. This product is not considered a RCRA hazardous waste.
14. Transport information	
ADR/RID:	No
IMO/IMDG:	No
KAO/IATA:	No
15. Regulatory information	
SARA Section 355:	None of the ingredients is listed
SARA Section 313:	None of the ingredients is listed
California Proposition 65:	None of the ingredients is listed
16. Other information	
Puritan Medical Products Company LLC provides the information in this document in good faith and believes the information to be accurate. The chemical, physical, and toxicological properties of this product have not been thoroughly investigated. It is the responsibility of the buyer to research and understand safe methods of handling, storing, and disposal of this product. Puritan Medical makes no warranty with respect to such information and assumes no liability for any loss or injury, which may result from the use of this information. It is the buyers responsibility to comply with local, state, and federal regulations concerning use and disposal of this product.	