



Puritan®
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

PurFlock Ultra® Flocked Swab - Dry Transport System

Standard tip, polystyrene handle



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. PurFlock Ultra® flocked swabs provide excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

Tip

- 0.200" Dia. (5.08 mm)
- 0.69" Length (17.526 mm)

Handle

- 0.100" Dia. (2.54 mm) tapered to 0.066" (1.676 mm)
- 6.01" Length (152.654 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Packaged in a reclosable dry transport tube for ease in transport to lab
- Double swab to accommodate multiple sampling sites
- Made in USA

Applications

- For collection and release of biological specimens
- Buccal cell collection
- Specimen collection
- Forensic evidence collection
- Rapid diagnostic collection
- DNA Testing



Puritan®

Ref #	Packaging
25-3306-2U BT Sterile	10/50/2 (10 boxes per case of 50 individually wrapped tubes per box, 2 swabs per tube)

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA
ISO 9001 2008 ISO 13485:2003 CE

Tel 800 321-2313 (US & Canada Only)
or 207-876-3311

Web puritanmedproducts.com
Email sales@puritanmedproducts.com