



Calgiswab® Calcium Alginate Tipped Applicator

Mini-tip Urethro-Genital Swab, Flexible Aluminum handle



Description

This specialty applicator is tipped with calcium alginate fiber, an organic material derived from seaweed and carded in Puritan's facility to assure conformance with medical device specifications. Diagnostic kit component.

Dimensions

5.437" Overall Length (138.1 mm)

Tip

- 0.078" Dia. (1.981 mm)
- 0.437" Length (11.1 mm)

Handle

- 0.025" Dia. (0.635 mm)
- 6.5" Length (165.1 mm)

Qualities

- Medical grade quality
- Calcium alginate tip contains no fatty acids, which may be inhibitory to many strains of bacteria
- Thin flexible aluminum handle for minimal patient discomfort
- Made in USA

Applications

- Flexible urethro-genital swab
- Specimen collection for testing and screening
- STD testing
- Pediatrics
- Ophthalmology
- Veterinary care



Puritan®

Ref #	Packaging
25-801 A 50 Sterile	10/100 (10 boxes of 50 individually wrapped packages, 1 p/pkg)