



# Calgiswab® Calcium Alginate Tipped Applicator

Regular tip, polystyrene handle



- Calcium alginate tip

## Description

This specialty applicator is tipped with calcium alginate fiber, an organic material derived from seaweed and carded in Puritan's facility to assure conformance with medical device specifications. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

## Dimensions

**6" Overall Length (152.4 mm)**

### Tip

- 0.203" Dia. (5.563 mm)
- 0.687" Length (17.45 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (147.925 mm)

## Qualities

- Medical grade quality
- Calcium alginate tip contains no fatty acids, which may be inhibitory to many strains of bacteria
- Made in USA

## Applications

- Specimen collection for testing and screening
- Wound care

Ref #	Packaging
25-806 1PA Sterile	10/100/1 (10 boxes per case, 100 individually wrapped packages per box, 1 p/pkg)



**Puritan®**