




# Cap-Shure® HydraFlock® Tipped Applicator

Elongated tip, Polystyrene handle



## Description

This absorbent flocked swab has an elongated tip and is bonded securely to a polystyrene handle by an aqueous based adhesive. Aerated tip protector designed for convenient, reliable specimen, DNA or evidence collection. Comes as a unit with a reclosable, ventilated cap for air drying and transporting samples.

## Dimensions

5.93" Overall Length (150.62 mm)

### Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

### Handle

- 0.100" Dia. (2.54 mm)
- 5.93" Length (147.625 mm)

## Qualities

- Medical grade quality
- Elongated flock tip
- Tip suspended away from sides
- Can lay flat to air dry
- Made in USA

## Applications

Diagnostic / Medical / Forensics & DNA

- Specimen collection
- DNA collection
- Evidence Collection

Ref #	Packaging
25-3206-H EC Sterile	10/50/1 (10 boxes per case, 50 individually wrapped packages per box, 1 p/pkg)



**Puritan®**