



Cap-Shure® Cotton Tipped Applicator

Regular tip, wood handle



Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. Aerated tip protector designed for convenient, reliable specimen, DNA or evidence collection. Comes as a unit with a reclosable, ventilated cap for air drying and transporting samples.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.188" Dia. (4.78 mm)
- 0.625" Length (15.875 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.812" Length (147.625 mm)

Qualities

- Medical grade quality
- Tip suspended away from sides
- Can lay flat to air dry
- Made in USA

Applications

Diagnostic / Medical / Forensics & DNA

- Specimen collection
- DNA collection
- Evidence Collection

Ref #	Packaging
25-806 1WC EC Sterile	10/50/1 (10 boxes per case, 50 individually wrapped packages per box, 1 p/pkg)



Puritan®