



Cotton Tipped Applicator

Aluminum handle

800 C

25-800 C 50



- Flexible aluminum handle

Description

This specialty applicator is tipped with high quality USP-grade cotton that is carded and coiled on site to maintain consistency. Thin flexible aluminum shaft for minimal patient discomfort. Available sterile or non-sterile.

Dimensions

	Overall Length	Tip (length x dia.)	Handle (length x dia.)
800 C	6.547" 165.1 mm	0.687 x 0.078" 17.45 x 1.981 mm	6.5 x 0.035" 165.1 x 1.981 mm
25-800 C 50	5.5" 139.7 mm		5.375 x 0.035" 136.525 x 0.889 mm

Qualities

- Medical grade USP cotton tip
- Small tip optimal for uses when a traditional size tip is too large
- Flexible aluminum handle
- Made in USA

Applications

- Forensic evidence collection
- Diagnostic specimen collection for testing and screening
- Pediatric ENT
- Eye microbiology
- Veterinary care

Ref #	Packaging
800 C Non-sterile	10/100 (10 bags of 100 swabs per bag)
25-800 C 50 Sterile	10/50/1 (10 boxes of 50 individually wrapped packages, 1 p/pkg)



Puritan®