




# Cotton Tipped Applicator

Polystyrene handle



## Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive. Available either non-sterile or sterile.

## Dimensions

**3" Overall Length (76.2 mm)**

### Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.449 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 2.937" Length (74.599 mm)

## Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- All purpose general swab
- Wound care
- Application of topical dressings
- Specimen collection
- Cell collection for DNA testing
- Cosmetics and hotel amenities
- Veterinary care

Ref #	Packaging
803-PC Non-Sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-803 2PC Sterile	10/100/2 (100 boxes of 100 individually wrapped packages, 2 p/pkg)



**Puritan®**