



PurFlock Ultra® Flocked Swab - Dry Transport System

Standard tip, polystyrene handle with breakpoint at 80 mm



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. PurFlock Ultra® flocked swabs provide excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

Tip

- 0.200" Dia. (5.08 mm)
- 0.69" Length (17.526 mm)

Handle

- 0.100" Dia. (2.54 mm) tapered to 0.066" (1.676 mm)
- 6.01" Length (152.654 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Packaged in a reclosable dry transport tube for ease in transport to lab
- Made in USA

Applications

- For collection and release of biological specimens
- Buccal cell collection
- Specimen collection
- Forensic evidence collection
- Rapid diagnostic collection
- DNA Testing

Ref #	Packaging
25-3306-U BT Sterile	5/100/1 (5 boxes, 100 tubes per box, 1 swab per tube)



Puritan®