



Cotton Tipped Applicator

Wood handle



Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch shaft by an aqueous based adhesive. Available sterile or non-sterile.

Dimensions

3" Overall Length (76.2 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 2.906" Length (73.812 mm)

Qualities

- Medical grade quality
- Standard size cotton tip
- Absorbent tip
- Soft smooth tip
- Made in USA

Applications

- General purpose swab
- Medication applicator
- Wound care
- First-aid kit component
- Application of topical dressings

Ref #	Packaging
803-WC Non-Sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-803 2WC Sterile	10/100/2 (10 boxes of 100 individually wrapped packages, 2 p/pkg)



Puritan®