



Cotton Tipped Applicator - Economy Grade

Standard Tip, Wood Handle



- Regular size cotton tip

- Biodegradable wood handle

Description

This economy grade cotton tipped applicator has a standard cotton tip and wood handle. The tip is comprised of imported economy medical grade cotton fiber. It's bonded to a naturally biodegradable wood handle and is perfect for medical, industrial, health and beauty environments for tasks like general cleaning, wound care, patient care, first aid kits, removing waxes, particulates, and residues from machinery, and applying salves.

Dimensions

3" Overall Length (76.2 mm)

Tip

- 0.188" Dia. (4.78mm)
- 0.625" Length (15.875 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 2.906" Length (73.812 mm)

Qualities

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges

Applications

- Patient care
- Wound care
- First aid kit component
- Health and beauty - waxing or applying salves
- Removing particulates and residues from machinery or surfaces
- General cleaning

Ref #	Packaging
803-WC HOSPITAL Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-803 2WC HOSPITAL Sterile	10/100/2 (10 boxes per case of 100 individually wrapped packages per box, 2 p/pkg)



HOSPITAL®