




Cotton Tipped Applicator

Polystyrene handle



Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive. Available either non-sterile or sterile.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.449 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip
- Made in USA

Applications

- All purpose general swab
- Wound care
- Application of topical dressings
- Specimen collection

Ref #	Packaging
806-PC Non-Sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-806 1PC Sterile	10/100/1 (100 boxes of 100 individually wrapped packages, 1 p/pkg)
25-806 2PC Sterile	10/100/2 (100 boxes of 100 individually wrapped packages, 2 p/pkg)



Puritan®