



Cotton Tipped Applicator - Economy Grade

Standard tip, Wood handle



-Standard size cotton tip

- Biodegradable wood handle

Description

This economy grade cotton tipped applicator has a standard cotton tip and wood handle. The tip is comprised of imported economy medical grade cotton fiber. It's bonded to a naturally biodegradable wood handle and is perfect for medical, industrial, health and beauty environments for task's like general cleaning, wound care, patient care, first aid kits, removing waxes, particulates, and residues from machinery, and applying salves.

Dimensions

6" Overall Length (76.2 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 5.812" Length (147.624mm)

Qualities

- Biodegradable handle
- Birch has high tensile strength
- Splinter free edges

Applications

- Patient and wound care
- First aid kit component
- Health and beauty - waxing or applying salves
- Removing particulates and residues from machinery or surfaces
- General cleaning



Ref #	Packaging
806-WC HOSPITAL Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-806 1WC HOSPITAL Sterile	10/100/1 (10 boxes per case, 100 individually wrapped packages per box, 1 p/pkg)
25-806 2WC HOSPITAL Sterile	10/100/2 (10 boxes per case, 100 individually wrapped packages per box, 2 p/pkg)