



# Cotton Tipped Applicator

Polystyrene handle



## Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive. Colored handle for color coded swab uses.

## Dimensions

6" Overall Length (152.4 mm)

### Tip

- 0.203" Dia. (5.156mm)
- 0.687" Length (17.45 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

## Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip
- Colored handle for coded uses
- Regular size tip
- Made in USA

## Applications

- DNA Collection
- Paternity testing
- Diagnostic testing
- Specimen collection

Ref #	Packaging
25-806 2PC BLUE Sterile	10/100/2 (10 boxes of 100 individually wrapped packages, 2p/pkg)
25-806 2PC PINK Sterile	
25-806 2PC YELLOW Sterile	



**Puritan®**