




Cotton Tipped Applicator

Small tip, Wood handle



Description

This cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. The small cotton tip is ideal for hard-to-reach or small areas.

Dimensions

5.875" Overall Length (149.225 mm)

Tip

- 0.141" Dia. (3.581 mm)
- 0.437 Length (11.1 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 5.812" Length (147.625 mm)

Qualities

- Sturdy wood handle, naturally biodegradable
- Soft smooth tip
- Neutral wrapper
- Made in USA

Applications

Medical

- General purpose swab
- Applying adhesives or other liquids
- Used for patient care or wound cleaning
- Applying medicines during surgery
- Applying topical dressings or ointments
- Features a small tip for hard-to-reach areas
- Convenient 2 swabs per sterile package

Ref #	Packaging
25-868 2WCS NEUTRAL Sterile	10/100/2 (10 boxes per case, 100 packages per box, 2 swabs per package)



Puritan®