




# Cotton Tipped Applicator

Hollow polypropylene handle



## Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to an anti-static polypropylene handle by an aqueous based adhesive.

## Dimensions

3" Overall Length (76.2 mm)

### Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.449 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 2.937" Length (74.6 mm)

## Qualities

- Regular cotton tip
- Anti-static polypropylene handle
- Made in USA

## Applications

Medical/Critical Environments & Industrial

- Multi-purpose general applicator
- For cleaning and removing residues
- For use in circuit board manufacturing



**Puritan®**

Ref #	Packaging
803-PPC Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)