



# HydraFlock® Flocked Swab - Dry Transport System

Mini-tip, polystyrene handle



## Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. The patented HydraFlock® flocked swab provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

## Dimensions

### Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

### Tip

- 0.135" Dia. (3.429 mm)
- 0.800" Length (20.32 mm)

### Handle

- 0.098" Dia. (2.489 mm) tapered to 0.049" (1.245 mm)
- 5.812" Length (147.62 mm)

## Qualities

- Medical grade quality
- Highly absorbent capacity
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

## Applications

- For nasopharyngeal and urethral cell collection
- Collection of bacterial or microorganisms
- Forensic evidence collection
- Specimen collection

| Ref #                   | Packaging  |
|-------------------------|--|
| 25-3316-H BT<br>Sterile | 5/100/1<br>(5 boxes, 100 tubes per box, 1 swab per tube) |



**Puritan®**