



# Cotton Tipped Applicator

Wood handle, large tip



## Description

This extra-absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

## Dimensions

3" Overall Length (76.2 mm)

### Tip

- 0.234" Dia. (5.943 mm)
- 0.687" Length (17.449 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 2.906" Length (73.812 mm)

## Qualities

- Medical grade quality
- Extra absorbent, large tip
- Soft smooth tip
- Made in USA

## Applications

### Medical:

- Specimen collection
- General purpose swab
- Medication applicator
- Wound care
- Applying salves or ointments

### Critical Environment/Industrial:

- Applying adhesives or liquids
- Removing dirt, contamination, fluxes or resins
- Cleaning with IPA or other solvents

Ref #	Packaging
803-WCL Non-sterile	10/5/50 (10 boxes per case, 5 bags per box, 50 swabs per bag)



**Puritan®**