




# Cotton Tipped Applicator

Unbreakable plastic handle



## Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a unbreakable plastic handle by an aqueous based adhesive.

## Dimensions

4" Overall Length (101.6 mm)

### Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.449 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 3.875" Length (98.425 mm)

## Qualities

- Regular cotton tip
- Unbreakable plastic handle
- Made in USA

## Applications

Medical/Critical Environments & Industrial

- Multi-purpose general applicator
- Cleaning swab



**Puritan®**

Ref #	Packaging
804-PC UNBRK Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)