



Cotton Tipped Applicator

Polystyrene handle, large tip



Description

This extra-absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

Dimensions

6" Overall Length (76.2 mm)

Tip

- 0.25" Dia. (6.35 mm)
- 0.687" Length (17.449 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Extra absorbent, large tip
- Soft smooth tip
- Semi-flexible plastic handle
- Made in USA

Applications

- All purpose general swab
- Wound care
- Application of topical dressings
- Specimen collection

Ref #	Packaging
806-PCL Non-sterile	10/5/100 (10 boxes per case, 5 bags per box, 100 swabs per bag)



Puritan®