



# Cotton Tipped Applicator

Wood handle, large tip



## Description

This extra absorbent large cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch handle by an aqueous based adhesive.

## Dimensions

6" Overall Length (76.2 mm)

### Tip

- 0.234" Dia. (5.94 mm)
- 0.687" Length (17.449 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.624 mm)

## Qualities

- Medical grade quality
- Extra absorbent, large tip
- Soft smooth tip
- Made in USA

## Applications

Medical/Critical Envrionments & Industrial

- Applying salves or ointments
- Wound care
- Specimen collection
- Applying adhesives or other liquids
- Removing dirt, contamination, fluxes, or resins
- Cleaning with IPA or other solvents

Ref #	Packaging
806-WCL Non-sterile	10/5/100 (10 boxes per case, 5 bags per box, 100 swabs per bag)



**Puritan®**



# Cotton Tipped Applicator

Wood handle, large tip

Ref #: 806-WCL

## NVR

NON-VOLATILE RESIDUE

Test method: IDEMA method M7-98

NVR Solvent - Hexane

NVR 12 µg/swab

The minimum detection limit for NVR is 10 µg/swab

## ION CHROMATOGRAPHY TESTING RESULTS

IONIC CLEANLINESS

Test method: IDEMA M13-99 (anions) M12-99 (cations)

### Anions

Fluoride	ND
Chloride	1.1000 µg/swab
Nitrate	0.1080 µg/swab
Bromide	1.14000 µg/swab
Nitrate	0.2680 µg/swab
Sulfate	7.0000 µg/swab
Phosphate	5.0000 µg/swab

### Cations

Lithium	ND
Sodium	21.0000 µg/swab
Ammonium	1.8800 µg/swab
Potassium	8.0800 µg/swab
Magnesium	1.2800 µg/swab
Calcium	3.4400 µg/swab

ND - Non Detected

Test procedures available upon request.



**Puritan**