



Cotton Tipped Applicator

Tapered mini-tip, wood handle



Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. Tapered mini cotton tip ideal for hard-to-reach or small areas.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.188 Dia. (4.775 mm) Regular Tip
- 0.625" Length (15.88 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 5.812" Length (147.624 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Tapered mini-tip
- Made in USA

Applications

- Multi-purpose use
- Ideal for hard-to-reach areas
- General cleaning swab

Ref #	Packaging
821-WC Non-sterile	10/5/100 (10 boxes per case, 5 bags per box, 100 swabs per bag)



Puritan®