



Cotton Tipped Applicator

Double tipped (one end tapered), Wood handle



Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. It is double tipped, one end with a mini tapered cotton tip and the other with a regular cotton tip.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.085" Dia. (2.159 mm) Tapered Tip
- 0.188" Dia. (4.775 mm) Regular Tip
- 0.625" Length (15.88 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 5.812" Length (147.624 mm)

Qualities

- Medical grade quality
- Double ended
- One end tapered, one end regular
- Made in USA

Applications

Medical/ Critical Environments & Industrial

- General patient care
- General cleaning swab
- Removing dirt, contamination, fluxes or resins

Ref #	Packaging
821-WC DBL Non-Sterile	10/5/100 (10 boxes per case, 5 bags per box, 100 swabs per bag)



Puritan®