




Cotton Tipped Applicator

Tapered mini-tip, wood handle



Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. Tapered mini cotton tip ideal for hard-to-reach or small areas.

Dimensions

2.969" Overall Length (75.413 mm)

Tip

- 0.12" Dia. (1.27 mm)
- 0.5" Length (12.7 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 2.906" Length (73.812 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Tapered mini-tip
- Made in USA

Applications

- Multi-purpose use
- Ideal for hard-to-reach areas
- General cleaning swab

Ref #	Packaging
823-WC Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



Puritan®