




# Cotton Tipped Applicator

Double tapered mini-tip, wood handle



## Description

This cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy white birch handle by an aqueous based adhesive.

Double tapered mini cotton tip ideal for hard-to-reach or small areas.

## Dimensions

6" Overall Length (152.4 mm)

### Tip

- 0.120" Dia. (1.27 mm) Tapered to 0.050" (1.27mm)
- 0.5" Length (12.7 mm)

### Handle

- 0.083" Dia. (2.108 mm)
- 2.906" Length (73.812 mm)

## Qualities

- Medical grade quality
- Double tapered mini-tip
- Rigid wood handle
- Low lint tips
- Made in USA

## Applications

Critical Environments / Industrial / Medical

- Multi-purpose use
- Wound care
- General cleaning swab
- Removing epoxy residues from connectors
- Cleaning tubes, scopes, sensors, and other components

Ref #	Packaging
826-WC DBL Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



**Puritan®**