



Cotton Tipped Applicator

Small tip, wood handle



Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is packed in a vial pack, suitable for kits.

Dimensions

3" Overall Length (76.2 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 2.906" Length (73.812 mm)

Qualities

- Small cotton tip
- Rigid wood handle
- Vial pack
- Made in USA

Applications

Critical Environments / Industrial / Medical

- All purpose swab used for wound care and dressing
- Suitable for kits
- Used for general cleaning and cosmetics

Ref #	Packaging
833-WCS Non-sterile	50/100 (50 vials, 100 swabs per vial)



Puritan®