



Cotton Tipped Applicator

Small tip, wood handle



Description

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a naturally biodegradable wooden handle with a hot-melt glue adhesive.

Dimensions

3.032" Overall Length (77.013 mm)

Tip

- 0.156" Dia. (3.962 mm)
- 0.625" Length (15.875 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 2.906" Length (73.812 mm)

Qualities

- Small cotton tip
- Rigid wood handle
- Made in USA

Applications

Critical Environments / Industrial / Medical

- All purpose swab used for general cleaning of small areas
- Wound care and dressing
- Applying medicines

Ref #	Packaging
838-WCS Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



Puritan[®]