




Cotton Tipped Applicator

Double ended standard tip, Polystyrene handle



Description

This absorbent double cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to an semi-flexible plastic handle by an aqueous based adhesive.

Dimensions

6.125" Overall Length (155.575 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Standard, double cotton tip
- Semi-flexible, hollow polypropylene handle with cotton tip at each end for twice the use
- Made in USA

Applications

Medical / Critical Environments / Industrial

- General purpose swab
- Fiber optics component manufacturing
- For use in rework and repair
- For use in removing excess solder and fluxes from circuit boards
- Cleaning metal finishes

Ref #	Packaging
861-PC DBL Non-sterile	10/5/100 (10 boxes per case, 5 bags per box, 100 swabs per bag)



Puritan®