




Cotton Tipped Applicator

Double ended standard tip, Wood handle



Description

This absorbent double cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a wood handle by an aqueous based adhesive.

Dimensions

6.062" Overall Length (153.975 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.625 mm)

Qualities

- Medical grade quality
- Standard, double cotton tip
- Sturdy wood handle
- Made in USA

Applications

Medical / Critical Environments / Industrial

- General purpose swab
- Fiber optics component manufacturing
- For use in rework and repair
- For use in removing excess solder and fluxes from circuit boards
- Cleaning metal finishes

Ref #	Packaging
861-WC DBL Non-sterile	10/5/100 (10 boxes per case, 5 bags per box, 100 swabs per bag)



Puritan®