




# Cotton Tipped Applicator

Small tip, Polystyrene handle



## Description

This cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled on-site to maintain consistency. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive. The small cotton tip is ideal for hard-to-reach or small areas.

## Dimensions

6" Overall Length (152.4 mm)

### Tip

- 0.172" Dia. (4.269 mm)
- 0.687" Length (17.45 mm)

### Handle

- 0.100" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

## Qualities

- Semi-flexible, hollow polystyrene handle
- Features minimal cotton for small or hard-to-reach areas
- Made in USA

## Applications

Critical Environments / Industrial / Medical

- All purpose swab used for general cleaning of small areas
- Removing oil or dirt from machine parts

Ref #	Packaging
868-PCS Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



**Puritan®**