




Cotton Tipped Applicator

Small tip, Wood handle



Description

This cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled on-site to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. The small cotton tip is ideal for hard-to-reach or small areas.

Dimensions

5.875" Overall Length (149.225 mm)

Tip

- 0.141" Dia. (3.581 mm)
- 0.437" Length (11.1 mm)

Handle

- 0.083" Dia. (2.108mm)
- 5.812" Length (147.625 mm)

Qualities

- Sturdy wood handle
- Minimal cotton tip for small or hard-to-reach areas
- Made in USA

Applications

Critical Environments / Industrial / Medical / Veterinary

- All purpose swab used for general cleaning of small areas
- Removing oil or dirt from machine parts
- Cleaning sensors and components
- Cleaning medical devices
- Veterinary care

Ref #	Packaging
868-WCS Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



Puritan®