



HydraFlock® Flocked Swab - Dry Transport System

Mini-tip, polystyrene handle with breakpoint at 80 mm



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. The patented HydraFlock® flocked swab provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

Tip

- 0.135" Dia. (3.429 mm)
- 0.800" Length (20.32 mm)

Handle

- 0.098" Dia. (2.489 mm) tapered to 0.049" (1.245 mm)
- 6.0" Length (152.4 mm)

Qualities

- Handle with breakpoint at 80 mm
- Medical grade quality
- Highly absorbent capacity
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Applications

- For nasopharyngeal and urethral cell collection
- Collection of bacteria or microorganisms
- Forensic evidence collection
- Specimen collection

Ref #	Packaging
25-3316-H EMB 80 BT Sterile	5/100/1 (5 boxes, 100 tubes per box, 1 swab per tube)



Puritan®