




# Histobrush® Cytology Collection Device

Tapered tip, polystyrene handle



-Soft, nylon bristles

- Polystyrene handle

\*Image not actual size

## Description

The Puritan Histobrush® is made with soft, tapered nylon bristles that is gentle on patients and reduces discomfort. The soft nylon bristles are used for cervical specimen collection and pap smear tests.

## Dimensions

8" Overall Length (203.2 mm)

### Tip

- 0.200" Dia. (5.08 mm)
- 0.750" Length (19.05 mm)

### Handle

- 0.117" Dia. (2.971 mm)
- 7.219" Length (183.362 mm)

## Qualities

- Sturdy polystyrene handle
- Soft nylon bristles
- Tapered tip

## Applications

- Cervical specimen collection
- Pap smear tests
- Buccal cell collection
- Veterinary

Ref #	Packaging
2188 Non-sterile	10/100 (10 boxes per case, 100 per box)
25-2188 Sterile	10/100/1 (10 boxes per case, 100 individually wrapped packages per box, 1 p/pkg)



**Puritan®**