



Rovers® Cervex-Brush®

Polyethylene handle



- Polyethylene bristles, various lengths

- Polyethylene handle

*Image not actual size

Description

Cell collection with the Cervex-Brush® provides a clear and accurate sample which leads to early detection of cervical cancer. The reliability of smears obtained by the Cervex-Brush® reduces the need for repeat testing. The brush is made up of various lengths of semi-circular plastic bristles. The shape is such that the top edges of the brush follow the contours of the cervix.

Dimensions

8" Overall Length (203.2 mm)

Qualities

- Sturdy polyethylene handle
- Soft flexible brush tip
- All necessary cells can be collected in one movement
- Made in USA

Applications

- For cell collection to detect cervical cancer
- Pap smear specimen collection

Ref #	Packaging
2195 Non-sterile	35/25 (35 bags per case, 25 per bag)



Puritan®