




CerviSoft® Cytology Collection Device

Foam tip, polystyrene handle



-Non-abrasive foam tip

- Polystyrene handle

*Image not actual size

Description

The Puritan CerviSoft® cytology collection device is made with a non-abrasive, soft CerviSoft® tip that is gentle on patients and reduces discomfort. The foam head readily releases cervical cells onto slides for testing. Ideal for use in sexual assault kits.

Dimensions

7.062" Overall Length (179.375 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 1" Length (25.4 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 7.010" Length (178.054 mm)

Qualities

- Sturdy polystyrene handle
- Soft, non-abrasive foam tip
- Foam tip readily releases specimens onto slides

Applications

- Cervical specimen collection
- Sexual assault kit component
- DNA collection

Ref #	Packaging
2197 Non-sterile	10/100 (10 boxes per case, 100 per box)
25-2197 Sterile	10/100/1 (10 boxes per case, 100 individually wrapped packages per box, 1 p/pkg)



Puritan®