



Rovers® Cervex-Brush® Combi

Polyethylene handle



Description

The Cervex-Brush® Combi will collect both ectocervical and endocervical cells from patients. The convenient Cervex-Brush® Combi provides health professionals with the option of a simultaneous Pap smear collection with a single device insertion.

Dimensions

8" Overall Length (203.2 mm)

Qualities

- Sturdy polyethylene handle
- Soft flexible bristles do not collect wet, humid cells allowing for easy transfer of cytological material to a glass slide or preservative fluid
- Detachable tip allows for use with both liquid based and conventional cervical cytology tests
- Requires only 2 full clockwise rotations
- High cellular yield

Applications

- For cell collection of both ectocervical and endocervical cells for Pap smears

Ref #	Packaging
2195 COMBI Non-sterile	4/25 (4 bags per case, 25 per bag)



Puritan®