



Puritan®
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

FAB-SWAB PurFlock Ultra® Tipped Applicator

Elongated tip, polystyrene handle



Description

This DNA controlled PurFlock Ultra® flocked swab has a greater ability to collect, retain and elute specimens than traditional spun swabs. The convenient, reclosable, dry transport tube is designed for reliable sample and evidence collection. The patented tube design entails a breathable filter label that covers holes to prohibit mold growth during transport. ID labels are included.

All of our DNA controlled swabs are manufactured in a DNA controlled and monitored clean room environment.

Dimensions

6" Overall Length (151.76 mm)

Tip

- 0.219" Dia. (5.56 mm)
- 0.685" Length (17.40 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.067" (1.702 mm)
- 5.96" Length (151.38 mm)

Qualities

- Medical grade quality
- Produced in a DNA controlled environment
- Reclosable transport tube
- Elongated flock tip
- Can be used dry, with IPA or sterile water
- Includes evidence information label
- Individually wrapped
- Made in USA

Applications

- DNA testing
- Forensic evidence collection



Puritan®

Ref #	Packaging
25-3606-U BT FABUSA Sterile	10/50/1 (10 boxes per case, 50 individually wrapped tubes per box, 1 swab per tube)

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA
ISO 9001 2008 ISO 13485:2003 CE

Tel 800 321-2313 (US & Canada Only)
or 207-876-3311

Web puritanmedproducts.com
Email sales@puritanmedproducts.com