





FAB-SWAB Cotton Tipped Applicator

Regular tip, polystyrene handle



Description

This DNA controlled absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a semi-flexible polystyrene handle.

The convenient, reclosable, dry transport device is designed for reliable sample and evidence collection. The patented tube design entails a breathable filter label that covers holes to prohibit mold growth during transport. ID labels are included.

Dimensions

5" Overall Length (127 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 4.875" Length (123.825 mm)

Qualities

- Medical grade quality
- Produced in a DNA controlled environment
- Reclosable transport tube
- Can be used dry, with IPA or sterile water
- Includes evidence information label
- Individually wrapped
- Made in USA

Applications

- DNA testing
- Forensic evidence collection

Ref #	Packaging
25-805 1PC TT FABUSA Sterile	10/50/1 (10 boxes per case, 50 individually wrapped tubes per box, 1 swab per tube)



Puritan®