



# FAB-SWAB Cotton Tipped Applicator

Regular tip, polystyrene handle



## Description

This DNA controlled absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a semi-flexible polystyrene handle.

The convenient, reclosable, dry transport device is designed for reliable sample and evidence collection. The patented tube design entails a breathable filter label that covers holes to prohibit mold growth during transport. ID labels are included.

## Dimensions

5.937" Overall Length (150.88 mm)

### Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

### Handle

- 0.100" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

## Qualities

- Medical grade quality
- Produced in a DNA controlled environment
- Reclosable transport tube
- Can be used dry, with IPA or sterile water
- Includes evidence information label
- Individually wrapped
- Made in USA

## Applications

- DNA testing
- Forensic evidence collection

Ref #	Packaging
25-806 1PC BT FABUSA Sterile	10/50/1 (10 boxes per case, 50 individually wrapped tubes per box, 1 swab per tube)



**Puritan®**