




FAB-MINI-AP Cotton Tipped Applicator

Mini-tip tip, wood handle



Description

This DNA controlled absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a wood handle.

The convenient, reclosable, dry transport device is designed for reliable sample and evidence collection. The patented tube design entails a breathable filter label that covers holes to prohibit mold growth during transport. ID labels are included.

Dimensions

5.095" Overall Length (129.413 mm)

Tip

- 0.12" Dia. (3.048 mm) tapered to 0.05" (1.27 mm)
- 0.5" Length (12.7 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 5.032" Length (127.813 mm)

Qualities

- Medical grade quality
- Produced in a DNA controlled environment
- Reclosable transport tube
- Can be used dry, with IPA or sterile water
- Includes evidence information label
- Tapered mini-tip for reaching into small areas
- Individually wrapped
- Made in USA

Applications

- DNA testing
- Forensic evidence collection

| Ref # | Packaging |
|----------------------------------|--|
| 25-825 1WCS TT FABUSA Sterile | 10/50/1 (10 boxes per case, 50 individually wrapped tubes per box, 1 swab per tube) |



Puritan®