




DNA-Free Foam Tipped Applicator

Large round tip, polypropylene handle



Description

This foam tipped applicator is produced with 100ppi medical grade foam. It has a large round head that is excellent for maximum cell collection.

Tested by an independent lab to the most stringent standards (less than 25 pg of human DNA), certified to be DNA-free, DNase-free and RNase-free

Dimensions

5.125" Overall Length (130.18 mm)

Tip

- 0.906" Dia. (23.0 mm)
- 1.125" Length (28.58 mm)
- 0.250" Thick (6.35 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5" Length (127.0 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip
- Made in USA

Applications

- Applying salves or ointments
- Wound care
- Application of topical dressings
- Specimen collection
- Applying adhesives or other liquids
- Removing dirt, contamination, fluxes or resins



Ref #	Packaging
25-1805 1PF RND FDNA Sterile	10/50/1 (10 boxes per case, 50 packages per box, 1 swab per package)