




DNA-Free PurFlock Ultra Tipped Applicator

Standard tip, polystyrene handle with molded breakpoint



Description

PurFlock's purity make it ideal for PCR, molecular assays, DNA testing, forensic applications as well as direct antigen testing. Its uniquely designed fiber structure accounts for better cell yield in the specimen where sample size is limited.

Tested by an independent lab to the most stringent standards (less than 25 pg of human DNA), certified to be DNA-free, DNase-free and RNase-free

Dimensions

6.04" Overall Length (153.416 mm)

Tip

- 0.200" Dia. (5.08 mm)
- 0.69" Length (17.53 mm)

Handle

- 0.1" Dia. tapered to 0.066" (2.54 mm to 1.68 mm)
- 6.01" Length (152.65 mm)

Qualities

- Medical grade quality
- Highly absorbent flock tip
- Made in USA

Applications

- Forensic / DNA testing
- Crime scene investigating
- Cell collection for DNA testing
- Specimen collection



Ref #	Packaging
25-3306-U TT FDNA Sterile	10/50/1 (10 boxes per case, 50 individually wrapped tubes per package, 1 swab per tube)